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| ANTIGEN           | BIOLOGICAL IMPORTANCE  | DIAGNOSTIC RELEVANCE  | MIKROGEN FEATURE   |
|-------------------|--|---|--|
| VlsE              | <ul style="list-style-type: none"> <li>• Variable membrane protein containing highly conserved regions</li> <li>• Circumvention of the host immune system</li> </ul> | <ul style="list-style-type: none"> <li>• Highly specific and sensitive</li> <li>• Key antigen for IgG detection</li> <li>• Early marker</li> </ul>  | <ul style="list-style-type: none"> <li>• Recombinant</li> <li>• Fusion protein from several genospecies</li> </ul> |
| OspC              | <ul style="list-style-type: none"> <li>• Surface protein</li> <li>• Important for the passage from tick to human</li> <li>• Binds plasminogen</li> </ul>             | <ul style="list-style-type: none"> <li>• Highly specific and sensitive</li> <li>• Key antigen for IgM detection</li> <li>• Early marker</li> </ul>  | <ul style="list-style-type: none"> <li>• Recombinant</li> <li>• 4 genospecies</li> </ul>                           |
| p18 (DbpA, Osp17) | <ul style="list-style-type: none"> <li>• Decorin-binding protein</li> <li>• Essential for spreading of Borrelia around the body</li> </ul>                           | <ul style="list-style-type: none"> <li>• Most heterogenous antigen between species</li> <li>• Highly specific and, when suitably combined, highly sensitive</li> <li>• Key antigen for IgG detection</li> </ul> | <ul style="list-style-type: none"> <li>• Recombinant</li> <li>• 5 genospecies</li> </ul>                           |
| p100              | <ul style="list-style-type: none"> <li>• Membrane protein</li> <li>• Function unknown</li> </ul>   | <ul style="list-style-type: none"> <li>• Highly specific</li> <li>• IgG marker</li> </ul>   | <ul style="list-style-type: none"> <li>• Recombinant</li> </ul>  |
| OspA              | <ul style="list-style-type: none"> <li>• Surface protein</li> <li>• Binds to the intestinal tract of the tick</li> </ul>   | <ul style="list-style-type: none"> <li>• (Detection of immunised people (USA))</li> </ul>   | <ul style="list-style-type: none"> <li>• Recombinant</li> </ul>  |
| p39 (BmpA)        | <ul style="list-style-type: none"> <li>• Borrelial membrane protein</li> </ul>   | <ul style="list-style-type: none"> <li>• Very specific and sensitive IgG marker</li> </ul>  | <ul style="list-style-type: none"> <li>• Recombinant</li> </ul>  |
| p58 (OppA-2)      | <ul style="list-style-type: none"> <li>• Oligopeptide-binding protein</li> </ul>   | <ul style="list-style-type: none"> <li>• Very specific and sensitive IgG marker</li> <li>• Scientifically evaluated</li> </ul>  | <ul style="list-style-type: none"> <li>• Recombinant</li> </ul>  |
| p41 (Fla)         | <ul style="list-style-type: none"> <li>• Component of flagellin</li> </ul>   | <ul style="list-style-type: none"> <li>• Relatively unspecific</li> </ul>   | <ul style="list-style-type: none"> <li>• Recombinant</li> </ul>  |

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Source : MIKROGEN GmbH, Neuried, Germany